




CoFlex[®]·TLC Lite

2 LAYER COMPRESSION

Provides Lite Compression
Patients with ABI ≥ 0.5

CoFlex TLC Lite 2 Layer Compression

A 2 step compression system designed to treat venous leg ulcers and other related conditions.



Both layers are **not made with natural rubber latex** and bond firmly together when applied which keeps the system in place.

(Patent Pending)

ORDERING INFORMATION		
Catalog # 7802TLC-TN	Rolls/Box 2	Boxes/Case 8
<i>Each Box Contains 2 Rolls - One of Each Layer</i>		
Roll	Description	Roll Size
Layer 1	Absorptive Foam Layer	10 cm x 3,1 m (Unstretched)
Layer 2	Compression Layer	10 cm x 4,7 m (Stretched)

Nylon Stocking included in each kit to apply over the completed dressing for patient comfort and ease of movement under clothes and on bed sheets.

Distributed By:

Sourcemed
surgicals | wound care

+27 (0)21 461 2985 www.sourcemed.co.za

